

Biotin Rat anti Mouse CD185/CXCR5 mAb

Catalog No.: A27444

Basic Information

Observed MW

Calculated MW

42kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0042

Conjugate

Biotin

Recommended Dilutions

FC 1 µg per 10⁶ cells in
100 µl volume

Background

Predicted to enable C-C chemokine binding activity and C-C chemokine receptor activity. Acts upstream of or within leukocyte chemotaxis and lymph node development. Located in external side of plasma membrane. Is expressed in several structures, including central nervous system; extraembryonic component; liver; and retina. Orthologous to human CXCR5 (C-X-C motif chemokine receptor 5).

Immunogen Information

Gene ID

12145

Swiss Prot

Q04683

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Blr1; Gpcr6; MDR15; CXC-R5; CXCR-5

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

IgG2a

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, pH7.2.

Validation Data

