

APC/Cyanine7 Rat anti-Mouse Ly-6C mAb

Catalog No.: A27499

Basic Information

Observed MW

Calculated MW

14kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ATC0017

Conjugate

APC-Cy7. Ex:651nm. Em:779nm.

Background

Predicted to enable acetylcholine receptor binding activity and acetylcholine receptor inhibitor activity. Predicted to be involved in acetylcholine receptor signaling pathway. Located in external side of plasma membrane. Is expressed in embryo mesenchyme; gut; liver; metanephros; and nasal cavity. Orthologous to human LY6H (lymphocyte antigen 6 family member H).

Recommended Dilutions

FC 0.25 µg per 10⁶ cells in
100 µl volume.

Immunogen Information

Gene ID

17067

Swiss Prot

POCW02

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Ly6c; Ly-6C; Ly-6C1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG2c, kappa

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, pH7.2.

Validation Data

