

# APC Rat anti-Mouse Ly6G mAb

Catalog No.: A27518

## Basic Information

---

### Observed MW

### Calculated MW

14kDa

### Category

Primary antibody

### Applications

FC

### Cross-Reactivity

Mouse

### CloneNo number

ATC0030

### Conjugate

APC. Ex:650nm. Em:660nm.

## Background

---

Located in external side of plasma membrane.

## Recommended Dilutions

---

FC 5  $\mu$ l per  $10^6$  cells in  
100  $\mu$ l volume

## Immunogen Information

---

### Gene ID

546644

### Swiss Prot

P35461

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

Gr1; Gr-1; Ly-6G

## Contact

---

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Rat

### Isotype

Rat / IgG2a, kappa

### Purification

Affinity purification


### Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  C57BL/6 mouse bone marrow cells were surface-stained with ABflo® 450 Rabbit anti-Mouse CD11b mAb (A26229, 5  $\mu$ l/Test) and APC Rat IgG2a isotype control (5  $\mu$ l/Test, left) or APC Rat anti-Mouse Ly6G mAb (A27518, 5  $\mu$ l/Test, right).