

# midP0032 Rabbit pAb

Catalog No.: A27532

## Basic Information

---

**Observed MW**

No

**Calculated MW**

No

**Category**

Primary antibody

**Applications**

WB

**Cross-Reactivity**

Human

## Background

---

No

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

**Gene ID**

No

**Swiss Prot**

No

**Immunogen**

This information is considered to be commercially sensitive.

**Synonyms**

No

## Contact

---

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Rabbit

**Isotype**

No

**Purification**

Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.