

Tcl1 Rabbit pAb

Catalog No.: A27540

Basic Information

Observed MW

Refer to figures

Calculated MW

14kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Mouse

Background

Predicted to enable identical protein binding activity; protein kinase binding activity; and protein serine/threonine kinase activator activity. Involved in activation of protein kinase B activity. Acts upstream of or within negative regulation of gene expression; regulation of stem cell population maintenance; and somatic stem cell population maintenance. Located in cell cortex; ooplasm; and pronucleus. Is expressed in several structures, including early conceptus; hemolymphoid system gland; liver; and seminiferous cord. Orthologous to human TCL1A (TCL1 family AKT coactivator A).

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

21432

Swiss Prot

P56280

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Tcl1a

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 50% glycerol, pH7.3.