

APC Rabbit anti-Human CD81 mAb

Catalog No.: A27605

Basic Information

Observed MW

Calculated MW

26kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC74354

Conjugate

APC. Ex:650nm. Em:660nm.

Background

The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell surface glycoprotein that is known to complex with integrins. This protein appears to promote muscle cell fusion and support myotube maintenance. Also it may be involved in signal transduction. This gene is localized in the tumor-suppressor gene region and thus it is a candidate gene for malignancies. Two transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

975

Swiss Prot

P60033

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

S5.7; CVID6; TAPA1; TSPAN28

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.


Validation Data




Flow cytometry: 1×10^6 HepG2 cells (negative control, left) and Jurkat cells (right) were surface-stained with APC Rabbit anti-Human CD81 mAb (A27605,5 µl/Test, orange line) or APC Rabbit IgG isotype control (A24173,5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1×10^6 Jurkat cells were surface-stained with APC Rabbit IgG isotype control (A24173,5 µl/Test, left) or APC Rabbit anti-Human CD81 mAb (A27605,5 µl/Test, right)



Flow cytometry: 1×10^6 Human PBMC were surface-stained with APC Rabbit anti-Human CD81 mAb (A27605,5 µl/Test, orange line) or APC Rabbit IgG isotype control (A24173,5 µl/Test, blue line). Non-fluorescently stained cells were used as blank control (red line). Cells in the lymphocyte gate were used for analysis.



Flow cytometry: 1×10^6 Human PBMC were surface-stained with APC Rabbit IgG isotype control (A24173,5 µl/Test, left) or APC Rabbit anti-Human CD81 mAb (A27605,5 µl/Test, right). Cells in the lymphocyte gate were used for analysis.