PE anti-Human Bcl-2 mAb

Catalog No.: A27722



Basic Information

Observed MW

Calculated MW

26kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

AMC0840

Conjugate

PE. Ex:565nm. Em:574nm.

Background

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

FC

5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID 596 **Swiss Prot**

P10415

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Bcl-2; PPP1R50

Contact

| a | 400-999-6126 |
|---------------------|---------------------------|
| \bowtie | cn.market@abclonal.com.cn |
| $\overline{\alpha}$ | www.abclonal.com.cn |

Product Information

SourceIsotypePurificationMouseMouse IgG1Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.2.