

PE anti-Mouse CD80/B7-1 Antibody

Catalog No.: A27786

Basic Information

Observed MW

Calculated MW

34kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ACC0017

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Predicted to enable coreceptor activity. Acts upstream of or within T cell costimulation; cellular response to lipopolysaccharide; and positive regulation of alpha-beta T cell proliferation. Located in external side of plasma membrane. Is expressed in heart ventricle and skin. Human ortholog(s) of this gene implicated in type 1 diabetes mellitus. Orthologous to human CD80 (CD80 molecule).

Recommended Dilutions

FC ≤0.2 µg per million cells
in 100 µl volume

Immunogen Information

Gene ID

12519

Swiss Prot

Q00609

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

B71; Ly53; TSA1; Cd28l; Ly-53; MIC17

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Armenian Hamster

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide and BSA,pH7.2.