

# Biotin Rabbit anti-Mouse IgA (Fc) mAb

Catalog No.: A28356

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Primary antibody

**Applications**

FC (intra)

**Cross-Reactivity**

Mouse

**CloneNo number**

ARC71353

**Conjugate**

Biotin

## Recommended Dilutions

---

**FC (intra)** ≤0.25 µg per million cells  
in 100 µl volume

## Background

---

## Immunogen Information

---

**Gene ID**

/

**Swiss Prot**

**Immunogen**

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

**Synonyms**

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

## Validation Data

---



Flow cytometry:  $1 \times 10^6$  293T (Transfection) cells were intracellularly-stained with Biotin Rabbit IgG isotype control (A25626, 5  $\mu$ l/Test, left) or Biotin Rabbit anti-Mouse IgA (Fc) mAb (A28356, 0.25  $\mu$ g, right), followed by APC Streptavidin staining.