

APC Rat anti-Mouse CD172a/SIRP α mAb

Catalog No.: A28435

Basic Information

Observed MW

Calculated MW

56 kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

Conjugate

APC. Ex:650nm. Em:660nm.

Background

Enables protein phosphatase binding activity; protein phosphatase inhibitor activity; and protein phosphorylated amino acid binding activity. Involved in several processes, including negative regulation of cytokine production; negative regulation of nitric oxide biosynthetic process; and negative regulation of signal transduction. Acts upstream of or within several processes, including hematopoietic progenitor cell differentiation; phagocytosis; and positive regulation of phagocytosis. Located in plasma membrane. Is expressed in several structures, including adrenal gland; cartilage; central nervous system; genitourinary system; and retina. Orthologous to human SIRPA (signal regulatory protein alpha).

Recommended Dilutions

FC $\leq 0.5 \mu\text{g}$ per million cells
in 100 μl volume

Immunogen Information

Gene ID

19261

Swiss Prot

P97797

Immunogen

This information is considered to be commercially sensitive.

Synonyms

Bit; P84; SIRP; SHP-1; CD172a; Ptpns1; SHPS-1; Idd13.2

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rat

Isotype

Rat / IgG1, kappa

Purification

Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with BSA, pH7.2.