

PE/Cyanine7 Rat anti-Mouse CD117/c-Kit mAb

Catalog No.: A28479

Basic Information

Observed MW

Calculated MW
109kDa

Category
Primary antibody

Applications
FC

Cross-Reactivity
Mouse

Conjugate
PE-Cy7. Ex:565nm. Em:778nm.

Background

The c-Kit proto-oncogene is the cellular homolog of the transforming gene of a feline retrovirus (v-Kit). The c-kit protein includes characteristics of a protein kinase transmembrane receptor. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

FC 1 µg per million cells in
100 µl volume

Immunogen Information

Gene ID	Swiss Prot
16590	P05532

Immunogen
This information is considered to be commercially sensitive.

Synonyms
W; Bs; Fdc; Ssm; SCO1; SCO5; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3

Contact

 | 400-999-6126
 | cn.market@abclonal.com.cn
 | www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Rat	Rat / IgG2b, kappa	Affinity purification

Storage
Store at 2-8°C. Avoid freeze.
Buffer: PBS with 0.09% Sodium azide and BSA ,pH7.2.