ABclonal www.abclonal.com

ABflo® 647 Rabbit anti-Human CD82 mAb

Catalog No.: A28484

Basic Information

Observed MW

Calculated MW 27kDa/30kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC77453

Conjugate

ABflo® 647. Ex:648nm. Em:664nm.

Background

This metastasis suppressor gene product is a membrane glycoprotein that is a member of the transmembrane 4 superfamily. Expression of this gene has been shown to be downregulated in tumor progression of human cancers and can be activated by p53 through a consensus binding sequence in the promoter. Its expression and that of p53 are strongly correlated, and the loss of expression of these two proteins is associated with poor survival for prostate cancer patients. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Recommended Dilutions

FC

5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID 3732 Swiss Prot

P27701

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

R2; 4F9; C33; IA4; ST6; GR15; KAI1; SAR2; TSPAN27

Contact

2		400-999-6126
\bowtie	Τ	cn.market@abclonal.com.cn
•	Т	www.abclonal.com.cn

Product Information


SourceIsotypePurificationRabbitIgGAffinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10^6 Human PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD3 mAb (A27177,5 μ I/Test) and ABflo® 647 Rabbit IgG isotype control (A22070,5 μ I/Test,left) or ABflo® 647 Rabbit anti-Human CD82 mAb (A28484,5 μ I/Test,right). Cells in the lymphocyte gate were used for analysis.