


Validation Data



Flow cytometry: 1×10^6 Human PBMC (untreated, left) and Human PBMC (treated with 20 ng/ml IL-4 and 20 ng/ml GM-CSF for 7 days, and 20 ng/ml TNF- α for last two days, right) were surface-stained with ABflo® 488 Mouse anti-Human CD11c mAb (A27558, 5 μ l/Test) and APC Rabbit anti-Human CD83 mAb (A28804, 5 μ l/Test). Cells in the monocyte gate were used for analysis.