

PE Rabbit anti-Mouse CD3e mAb

Catalog No.: A28812

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC80921

Conjugate

PE. Ex:565nm. Em:574nm.

Background

CD3ε is a 20 kD transmembrane protein, also known as CD3 or T3. It is a member of the Ig superfamily, mainly expressed on T cells and NK-T cells, and expressed at different levels on thymocytes during the process of T cell differentiation. CD3ε forms TCR complexes by binding to CD3δ, γ and δ chains as well as TCR α/β or γ/δ chains. CD3 plays a key role in TCR signal transduction, T cell activation and antigen recognition by binding to peptide /MHC antigen complexes.

Recommended Dilutions

FC ≤0.25 µg per million cells in
100 µl volume

Immunogen Information

Gene ID

12501

Swiss Prot

P22646

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CD3;CD3epsilon;T3e

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse splenocytes were surface-stained with ABflo® 450 Rabbit anti-Mouse CD19 mAb (A26192, 5 µl/Test) and PE Rabbit IgG isotype control (A24172, 5 µl/Test, left) or PE Rabbit anti-Mouse CD3e mAb (A28812, 0.25 µg, right).