

Rabbit anti-Human BTLA mAb

Catalog No.: A28976

Basic Information

Observed MW

Calculated MW

33kDa

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human

CloneNo number

ARC5251

Background

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.

Recommended Dilutions

FC 1:100-1:500

Immunogen Information

Gene ID

151888

Swiss Prot

Q7Z6A9

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

BTLA1; CD272

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: 25mM HEPES, 150mM NaCl, 0.03% Proclin, 50% glycerol, pH8.0.