

FITC Rabbit anti-Human/Monkey CD64 mAb

Catalog No.: A29096

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

FC

Cross-Reactivity

Human, Cynomolgus

CloneNo number

ARC82983

Conjugate

FITC. Ex:491nm. Em:516nm.

Background

This gene encodes a protein that plays an important role in the immune response. This protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family members located on chromosome 1.

Recommended Dilutions

FC 5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID

Hu 2209/Cyno 102147198

Swiss Prot

P12314

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CD64; CD64A; FCG1; FCGR1; FCRI; FcgammaRI; IGFR1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data



Flow cytometry: 1×10^6 Cynomolgus PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD14 mAb (A25196,5 μ l/Test) and FITC Rabbit IgG isotype control (A25616,5 μ l/Test, left) or FITC Rabbit anti-Human/Monkey CD64 mAb (A29096,5 μ l/Test, right). Cells in the monocyte gate were used for analysis.

Flow cytometry: 1×10^6 Human PBMC were surface-stained with ABflo® 450 Rabbit anti-Human/Monkey CD14 mAb (A25196,5 μ l/Test) and FITC Rabbit IgG isotype control (A25616,5 μ l/Test, left) or FITC Rabbit anti-Human/Monkey CD64 mAb (A29096,5 μ l/Test, right). Cells in the monocyte gate were used for analysis.