

FOXD3 Rabbit pAb

Catalog No.: A2926

Basic Information

Observed MW

48kDa

Calculated MW

48kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse

Background

This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

27022

Swiss Prot

Q9UJU5

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 322-335AA of human FOXD3.

Synonyms

AIS1; HFH2; VAMAS2; Genesis; FOXD3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from HeLa cells, using FOXD3 Rabbit pAb (A2926).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.