

# FOXN2 Rabbit pAb

Catalog No.: A2933

## Basic Information

### Observed MW

### Calculated MW

47kDa

### Category

Primary antibody

### Applications

ELISA,IF/ICC

### Cross-Reactivity

Human

## Background

This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

## Recommended Dilutions

IF/ICC 1:50 - 1:200

## Immunogen Information

### Gene ID

3344

### Swiss Prot

P32314

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 172-431 of human FOXN2 (NP\_002149.2).

### Synonyms

HTLF; FOXN2

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of A549 cells using FOXN2 antibody (A2933).