

FOXO3A Rabbit mAb

Catalog No.: A29347 **Recombinant**

Basic Information

Observed MW

82-97 kDa

Calculated MW

71 kDa/48 kDa

Category

Primary antibody

Applications

WB, IP, IHC-P, ELISA, WB

Cross-Reactivity

Human, Mouse

CloneNo number

ARC3594

Background

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.

Recommended Dilutions

WB 1:20000 - 1:40000

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

IHC-P 1:6000-1:30000

ELISA Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on your
specific assay requirements.

WB 1:3000 - 1:10000

Immunogen Information

Gene ID

2309

Swiss Prot

O43524

Immunogen

This information is considered to be commercially sensitive.

Synonyms

FOXO2; AF6q21; FKHL1; FOXO3A; FKHL1P2

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C


Contact

☎ | 400-999-6126


✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn


Validation Data



Western blot analysis of lysates from PC-3 cells using FOXO3A Rabbit mAb (A29347) at 1:20000 dilution incubated overnight at 4°C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5 s.




Western blot analysis of lysates from Mouse pancreas using FOXO3A Rabbit mAb (A29347) at 1:5000 dilution incubated overnight at 4°C. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 5 s.




Immunoprecipitation of FOXO3A from 300 µg extracts of PC-3 cells was performed using 2 µg of FOXO3A Rabbit mAb (A29347). Rabbit Control IgG (AC005) was used to precipitate the Control IgG sample. IP samples were eluted with 1x Laemmli Buffer. The Input lane represents 10% of the total input. Western blot analysis of immunoprecipitates was conducted using FOXO3A Rabbit mAb (A29347) at a dilution of 1:50000.


Validation Data



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma tissue using FOXO3A Rabbit mAb (A29347) at a dilution of 1:12000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using FOXO3A Rabbit mAb (A29347) at a dilution of 1:12000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Mouse spleen tissue using FOXO3A Rabbit mAb (A29347) at a dilution of 1:12000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.