

# AKAP8L Rabbit pAb

Catalog No.: A3050

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

72kDa

### Category

Primary antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human

## Background

Enables histone deacetylase binding activity and lamin binding activity. Involved in several processes, including mitotic chromosome condensation; regulation of histone modification; and regulation of mRNA export from nucleus. Located in chromatin; cytoplasm; and nuclear lumen. Part of ribonucleoprotein complex.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

26993

### Swiss Prot

Q9ULX6

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

HA95; HAP95; NAKAP; NAKAP95; AKAP8L

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.