

NARFL Rabbit pAb

Catalog No.: A3688 **1 Publications**

Basic Information

Observed MW

56kDa

Calculated MW

53kDa

Category

Primary antibody

Applications

WB,FC

Cross-Reactivity

Human

Background

Predicted to enable 4 iron, 4 sulfur cluster binding activity. Involved in several processes, including iron-sulfur cluster assembly; oxygen homeostasis; and response to hypoxia. Part of CIA complex.

Recommended Dilutions

WB	1:500 - 1:1000
FC	1:20 - 1:50

Immunogen Information

Gene ID

64428

Swiss Prot

Q9H6Q4

Immunogen

A synthetic peptide of human NARFL

Synonyms

PRN; HPRN; IOP1; NAR1; LET1L; NARFL

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

Validation Data



Western blot analysis of extracts of MDA-MB-435 cells, using NARFL antibody (A3688).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.