

# REEP5 Rabbit pAb

Catalog No.: A4191 **1 Publications**

## Basic Information

### Observed MW

20kDa

### Calculated MW

21kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human, Mouse

## Background

Predicted to be involved in endoplasmic reticulum organization and regulation of intracellular transport. Located in endoplasmic reticulum tubular network.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

7905

### Swiss Prot

Q00765

### Immunogen

A synthetic peptide of human REEP5

### Synonyms

DP1; TB2; YOP1; POB16; Yip2e; D5S346; C5orf18; REEP5

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of extracts of HeLa cells, using REEP5 antibody (A4191).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25 $\mu$ g per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.