

NGEF Rabbit pAb

Catalog No.: A4583

Basic Information

Observed MW

82kDa

Calculated MW

82kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Rat

Background

Enables guanyl-nucleotide exchange factor activity. Predicted to be involved in several processes, including activation of GTPase activity; ephrin receptor signaling pathway; and negative regulation of dendritic spine morphogenesis. Predicted to be located in cytosol. Predicted to be active in glutamatergic synapse.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

25791

Swiss Prot

Q8N5V2

Immunogen

A synthetic Peptide of human NGEF

Synonyms

EPHEXIN; ARHGEF27; NGEF

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.