

# HN1 Rabbit pAb

Catalog No.: A4716

## Basic Information

### Observed MW

28kDa

### Calculated MW

16kDa

### Category

Primary antibody

### Applications

WB,FC

### Cross-Reactivity

Human

## Background

Located in nuclear membrane; nucleolus; and nucleoplasm.

## Recommended Dilutions

**WB** 1:500 - 1:1000

**FC** 1:20 - 1:50

## Immunogen Information

### Gene ID

51155

### Swiss Prot

Q9UK76

### Immunogen

A synthetic Peptide of human HN1

### Synonyms

HN1; ARM2; HN1A

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.

## Validation Data

---



Western blot analysis of extracts of 293T cells, using HN1 antibody (A4716).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.