

IWS1 Rabbit pAb

Catalog No.: A4839

Basic Information

Observed MW

Calculated MW

92kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

Involved in regulation of histone modification; regulation of mRNA export from nucleus; and regulation of mRNA processing. Located in nucleoplasm.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55677

Swiss Prot

Q96ST2

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

IWS1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.