

FUT11 Rabbit pAb

Catalog No.: A5190

Basic Information

Observed MW

Calculated MW

56kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable alpha-(1->3)-fucosyltransferase activity. Predicted to be involved in fucosylation. Predicted to be located in Golgi cisterna membrane. Predicted to be integral component of membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

170384

Swiss Prot

Q495W5

Immunogen

Recombinant protein of human FUT11

Synonyms

FUCTXI; FUT11

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data

Western blot analysis of various lysates using FUT11 Rabbit pAb (A5190).

