

APOBEC3C Rabbit pAb

Catalog No.: A5877

Basic Information

Observed MW

Calculated MW

23kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

This gene is a member of the cytidine deaminase gene family. It is one of seven related genes or pseudogenes found in a cluster thought to result from gene duplication, on chromosome 22. Members of the cluster encode proteins that are structurally and functionally related to the C to U RNA-editing cytidine deaminase APOBEC1. It is thought that the proteins may be RNA editing enzymes and have roles in growth or cell cycle control.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

27350

Swiss Prot

Q9NRW3

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

A3C; PBI; ARP5; ARDC2; ARDC4; APOBEC1L; bK150C2.3; APOBEC3C

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.