

KDM3B/KDM3A Rabbit pAb

Catalog No.: A6158

Basic Information

Observed MW

Refer to figures

Calculated MW

84kDa/154kDa/191kDa/147kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
51780/55818

Swiss Prot
Q7LBC6/Q9Y4C1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms
KDM3B/KDM3A

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.