

UBE3C Rabbit pAb

Catalog No.: A6442 **2 Publications**

Basic Information

Observed MW

124kDa

Calculated MW

124kDa

Category

Primary antibody

Applications

ELISA, WB, IF/ICC

Cross-Reactivity

Human

Background

Enables ubiquitin protein ligase activity. Involved in protein polyubiquitination. Predicted to be located in nucleus. Predicted to be part of proteasome complex.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:10 - 1:100

Immunogen Information

Gene ID

9690

Swiss Prot

Q15386

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 970-1083 of human UBE3C (NP_055486.2).

Synonyms

RAUL; HECTH2; NDSMBA; NEDSMBA; UBE3C

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.