

POLM Rabbit pAb

Catalog No.: A6835

Basic Information

Observed MW

Refer to figures

Calculated MW

55kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable DNA-directed DNA polymerase activity. Predicted to be involved in double-strand break repair via nonhomologous end joining. Predicted to act upstream of or within B cell differentiation and somatic hypermutation of immunoglobulin genes. Predicted to be located in nucleoplasm. Predicted to be active in nucleus.

Recommended Dilutions

WB 1:200 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

27434

Swiss Prot

Q9NP87

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Tdt-N; Pol Mu; POLM

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.