

# ELF4 Rabbit pAb

Catalog No.: A7180

## Basic Information

### Observed MW

### Calculated MW

71kDa

### Category

Primary antibody

### Applications

ELISA,IF/ICC

### Cross-Reactivity

Human

## Background

The protein encoded by this gene is a transcriptional activator that binds and activates the promoters of the CSF2, IL3, IL8, and PRF1 genes. The encoded protein is involved in natural killer cell development and function, innate immunity, and induction of cell cycle arrest in naive CD8+ cells. Two transcript variants encoding the same protein have been found for this gene.

## Recommended Dilutions

IF/ICC 1:50 - 1:100

## Immunogen Information

### Gene ID

2000

### Swiss Prot

Q99607

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human ELF4 (NP\_001412.1).

### Synonyms

MEF; ELFR; AIFBL2; ELF4

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data

---



Immunofluorescence analysis of A549 cells using ELF4 antibody (A7180). Blue: DAPI for nuclear staining.