

MTMR4 Rabbit pAb

Catalog No.: A7761

Basic Information

Observed MW

Refer to figures

Calculated MW

133kDa

Category

Primary antibody

Applications

ELISA, IHC-P

Cross-Reactivity

Human, Mouse

Background

Enables protein phosphatase binding activity. Involved in regulation of phosphatidylinositol dephosphorylation. Located in endosome.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

9110

Swiss Prot

Q9NYA4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 570-820 of human MTMR4 (NP_004678.3).

Synonyms

ZFYVE11; FYVE-DSP2; MTMR4

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human hepatocholangiocarcinoma using MTMR4 Rabbit pAb (A7761) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded Human placenta using MTMR4 Rabbit pAb (A7761) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.