

SETD3 Rabbit pAb

Catalog No.: A8071 **2 Publications**

Basic Information

Observed MW

70kDa

Calculated MW

67kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

Enables protein-L-histidine N-tele-methyltransferase activity. Involved in actin modification and peptidyl-histidine methylation. Located in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:200

Immunogen Information

Gene ID

84193

Swiss Prot

Q86TU7

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SETD3 (NP_115609.2).

Synonyms

hSETD3; C14orf154; SETD3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.