

APBB1IP Rabbit pAb

Catalog No.: A8220

Basic Information

Observed MW

Refer to figures

Calculated MW

73kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to be involved in signal transduction. Predicted to act upstream of or within T cell activation via T cell receptor contact with antigen bound to MHC molecule on antigen presenting cell and positive regulation of cell adhesion. Located in cytosol and plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54518

Swiss Prot

Q7Z5R6

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RIAM; INAG1; PREL1; RARP1; APBB1IP

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.