

CCHCR1 Rabbit pAb

Catalog No.: A8221

Basic Information

Observed MW

Refer to figures

Calculated MW

89kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a protein with five coiled-coil alpha-helical rod domains that is thought to act as a regulator of mRNA metabolism through its interaction with mRNA-decapping protein 4. It localizes to P-bodies, the site of mRNA metabolism, with an N-terminus that is required for this subcellular localization, suggesting it is a P-body component. Naturally occurring mutations in this gene are associated with psoriasis.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

54535

Swiss Prot

Q8TD31

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HCR; SBP; pg8; C6orf18; CCHCR1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.