

HNRNPCL1 Rabbit pAb

Catalog No.: A8606

Basic Information

Observed MW

Refer to figures

Calculated MW

32kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human, Rat

Background

Enables identical protein binding activity. Predicted to be active in nucleus.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

343069

Swiss Prot

O60812

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HNRPCL1; HNRNPCL1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.