

# EPS8 Rabbit pAb

Catalog No.: A8826 **1 Publications**

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

92kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

This gene encodes a member of the EPS8 family. This protein contains one PH domain and one SH3 domain. It functions as part of the EGFR pathway, though its exact role has not been determined. Highly similar proteins in other organisms are involved in the transduction of signals from Ras to Rac and growth factor-mediated actin remodeling. Alternate transcriptional splice variants of this gene have been observed but have not been thoroughly characterized.

## Recommended Dilutions

WB 1:200 - 1:2000

## Immunogen Information

### Gene ID

2059

### Swiss Prot

Q12929

### Immunogen

A synthetic Peptide of human EPS8

### Synonyms

DFNB102; EPS8

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.