

# RAB27B Mouse mAb

Catalog No.: A8897 **1 Publications**

## Basic Information

### Observed MW

27kDa

### Calculated MW

25kDa

### Category

Primary antibody

### Applications

ELISA,WB,FC

### Cross-Reactivity

Human, Mouse

## Background

Enables guanyl ribonucleotide binding activity; myosin V binding activity; and protein domain specific binding activity. Involved in multivesicular body sorting pathway and positive regulation of exocytosis. Located in Golgi stack and cytoplasmic vesicle.

## Recommended Dilutions

WB	1:500 - 1:2000
FC	1:20 - 1:50

## Immunogen Information

### Gene ID

5874

### Swiss Prot

O00194

### Immunogen

A synthetic peptide of human RAB27B

### Synonyms

C25KG; RAB27B

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Mouse

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,pH7.3.