

# FLAD1 Rabbit pAb

Catalog No.: A8960

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

65kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

This gene encodes the enzyme that catalyzes adenylation of flavin mononucleotide (FMN) to form flavin adenine dinucleotide (FAD) coenzyme. Alternatively spliced transcript variants encoding distinct isoforms have been observed.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

80308

### Swiss Prot

Q8NFF5

### Immunogen

A synthetic Peptide of human FLAD1

### Synonyms

FAD1; FADS; PP591; LSMFLAD

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.