

TEKT4 Rabbit pAb

Catalog No.: A9241

Basic Information

Observed MW

Refer to figures

Calculated MW

51kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to be involved in cilium assembly and cilium movement involved in cell motility. Predicted to act upstream of or within regulation of brood size. Located in nucleus.

Recommended Dilutions

WB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

150483

Swiss Prot

Q8WW24

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TEKT4; tektin-4

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.