

USP28 Rabbit pAb

Catalog No.: A9292 **2 Publications**

Basic Information

Observed MW

150kDa

Calculated MW

122kDa

Category

Primary antibody

Applications

WB

Cross-Reactivity

Human

Background

The protein encoded by this gene is a deubiquitinase involved in the DNA damage pathway and DNA damage-induced apoptosis. Overexpression of this gene is seen in several cancers.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

57646

Swiss Prot

Q96RU2

Immunogen

A synthetic Peptide of human USP28

Synonyms

USP28

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of Jurkat cells, using USP28 antibody (A9292).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.