

# Rabbit Control IgG

Catalog No.: AC005N **1 Publications**

## Basic Information

---

**Observed MW**

**Calculated MW**

**Category**

Loading control antibody

**Applications**

IP

**Cross-Reactivity**

## Background

---

## Recommended Dilutions

---

IP 1:50 - 1:200

## Immunogen Information

---

**Gene ID** **Swiss Prot**

**Immunogen**

None

**Synonyms**

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Potein A/G purification


**Storage**

Store at 4°C Avoid freeze / thaw cycles.

Buffer: PBS,pH7.3.

## Validation Data

---



Immunoprecipitation of over-expressed DDDDK-tagged protein in 293T cells incubated using DDDDK-tag antibody (AE063). A mock served as negative control using rabbit Control IgG AC005N and over-expressed 293T cell lysate served as positive control.