

PE Mouse anti-HA.11 Epitope Tag mAb

Catalog No.: AE145

Basic Information

Observed MW

Calculated MW

Category

Tag antibody

Applications

FC

Cross-Reactivity

Species independent

CloneNo number

AMC0830

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Recommended Dilutions

FC 0.25 µg per 10⁶ cells in
100 µl volume

Immunogen Information

Gene ID

/

Swiss Prot

/

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HA epitope tag; HA1; HA2; hemagglutinin; Hemagglutinin HA1 chain; Hemagglutinin HA2 chain; YPYDVPDYA; Hemagglutinin tag

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Mouse

Isotype

Mouse IgG1, κ Isotype Ctrl

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,pH7.2.

Validation Data

