

Phospho-Bcl-2-S70 Rabbit pAb

Catalog No.: AP0315

Basic Information

Observed MW

Calculated MW

26kDa

Category

Primary antibody

Applications

Cross-Reactivity

Human

Background

This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

Immunogen Information

Gene ID

596

Swiss Prot

P10415

Immunogen

A phospho specific peptide corresponding to residues surrounding S70 of human Bcl-2

Synonyms

Bcl-2; PPP1R50; Phospho-Bcl-2-S70

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

