

Phospho-IRS1-S307 Rabbit pAb

Catalog No.: AP0371 **2 Publications**

Basic Information

Observed MW

Calculated MW

132kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a protein which is phosphorylated by insulin receptor tyrosine kinase. Mutations in this gene are associated with type II diabetes and susceptibility to insulin resistance.

Recommended Dilutions

Immunogen Information

Gene ID

3667

Swiss Prot

P35568

Immunogen

A phospho specific peptide corresponding to residues surrounding S307 of human IRS1

Synonyms

HIRS-1; Phospho-IRS1-S307

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.