

Phospho-IGF1R-Y1135/Y1136 + INSR-Y1150/Y1151 Rabbit pAb

Catalog No.: AP0547

Basic Information

Observed MW

Refer to figures

Calculated MW

154kDa/155kDa/156kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014]

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

3480/3643

Swiss Prot

P08069/P06213

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

IGF1R/INSR; Phospho-IGF1R-Y1135/Y1136 + INSR-Y1150/Y1151

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.