

# Phospho-RICTOR-T1135 Rabbit pAb

Catalog No.: AP0869

## Basic Information

### Observed MW

### Calculated MW

192kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

RICTOR and MTOR (FRAP1; MIM 601231) are components of a protein complex that integrates nutrient- and growth factor-derived signals to regulate cell growth (Sarbasov et al., 2004 [PubMed 15268862]).

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

253260

### Swiss Prot

Q6R327

### Immunogen

A synthetic phosphorylated peptide around T1135 of human RICTOR (NP\_689969.2).

### Synonyms

PIA; AVO3; hAVO3; Phospho-RICTOR-T1135

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.