

Phospho-TAZ-S89 Rabbit pAb

Catalog No.: AP0919 **2 Publications**

Basic Information

Observed MW

Refer to Figures

Calculated MW

44kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Enables transcription coactivator activity. Involved in several processes, including hippo signaling; positive regulation of cell differentiation; and regulation of signal transduction. Located in cytosol and nuclear body.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

25937

Swiss Prot

Q9GZV5

Immunogen

A synthetic phosphorylated peptide around S89 of human WWTR1 (NP_056287.1).

Synonyms

TAZ; Phospho-TAZ-S89

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.