

HRP-conjugated Rabbit anti-Sheep IgG (H+L)

Catalog No.: AS023 **3 Publications**

Basic Information

Observed MW

Calculated MW

Category

Secondary antibody

Applications

WB,IHC-P,ELISA

Cross-Reactivity

Conjugate

HRP

Background

Secondary antibodies are affinity-purified antibodies which will work with target-specific primary antibody in the detection, sorting or purification of its specified target. Secondary antibodies offer increased versatility enabling users to use many detection systems (e.g. HRP, AP, fluorescence). They can also provide greater sensitivity through signal amplification as multiple secondary antibodies. Most commonly, secondary antibodies are generated by immunizing the host animal (different from host species of primary antibody) with a pooled population of normal immunoglobulins from the host species of primary antibody and can be further purified and modified (i.e. antibody fragmentation, label conjugation, etc.) to ensure well-characterized specificity to corresponding normal immunoglobulins.

Recommended Dilutions

| | |
|--------------|-------------------|
| WB | 1:5000 - 1:100000 |
| IHC-P | 1:500 - 1:5000 |

Immunogen Information

| | |
|----------------|-------------------|
| Gene ID | Swiss Prot |
|----------------|-------------------|

Immunogen

Sheep IgG

Synonyms

Contact

| | |
|--|--------------|
|  | 400-999-6126 |
|--|--------------|

| | |
|--|---------------------------|
|  | cn.market@abclonal.com.cn |
|--|---------------------------|

| | |
|--|--|
|  | www.abclonal.com.cn |
|--|--|

Product Information

Source

Rabbit

Isotype

Horseradish peroxidase
conjugated IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.