ABclonal www.abclonal.com

ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb

Catalog No.: AS155

Basic Information

Background

Observed MW

Calculated MW

36kDa

Category

Secondary antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC68493

Conjugate

ABflo® 488. Ex:491nm. Em:516nm.

Recommended Dilutions

FC

5 μ l per 10^6 cells in 100 μ l volume

Immunogen Information

Gene ID

Swiss Prot P01868

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

a		400-999-6126
\bowtie		cn.market@abclonal.com.cn
\odot	Т	www.abclonal.com.cn

Product Information

SourceIsotypePurificationRabbitIgGAffinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide, 0.2% BSA, pH7.3.

Validation Data







Flow cytometry: 1X10^6 Raji cells (negative control,left) and Human PBMC (right) were surface-stained with Mouse anti Human CD3 mAb (2 µg/mL,orange line) or Mouse IgG isotype control (2 µg/mL,blue line), followed by ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb (AS155,5 µl/Test) staining. Nonfluorescently stained cells were used as blank control (red line).

Flow cytometry: 1X10^6 Human PBMC were surface-stained with Mouse IgG isotype control (2 μ g/mL,left) or Mouse anti Human CD3 mAb (2 μ g/mL,right). ABflo® 488 Rabbit anti-Mouse IgG (Fc) mAb (AS155,5 μ l/Test) was used as a secondary antibody.